



Transiti[©]ns[®] XTRActive[®]

A PHOTOCROMIC DYE CHEMIST'S CHALLENGE:

USING SCIENCE TO CREATE A PHOTOCROMIC
LENS THAT OFFERS EXTRA PROTECTION FROM
LIGHT OUTDOORS, INDOORS AND EVEN IN THE CAR



Access to
Quinn-Peterson

Access to
Quinn-Peterson


- Access to
Quinn-Peterson

Access to
Quinn-Peterson



A QUICK LOOK AT PHOTOCROMICS

(PHOTO = LIGHT, CHROMA = COLOR)



Photochromic technology has been a core competency at Transitions Optical for more than 25 years. The definition of photochromics is a light-induced, reversible change in color. A photochromic lens contains millions of photochromic molecules.

In the **colorless state**, these molecules consist of two smaller chromophores (or halves) held out of plane from one another (at different angles). Or, in other words, when the lenses are clear, the two halves of the photochromic molecule are perpendicular to each other.

When **exposed to light energy**, a chemical bond breaks and the molecule rearranges itself from two smaller chromophores into one large chromophore (or one flat plane) that absorbs in both the UV and the visible portion of the electromagnetic (EM) light spectrum.

This means the molecule now has color. Simply put, the UV light changes the shape of the photochromic molecules and the lenses darken. When the light **energy is taken away**, heat reforms the bond, and the lenses go back to clear.

Because the fade back reaction is driven by heat, when it is very hot outside the photochromic molecules are pushed back to their clear state even as UV light energy is driving them to activate. This is why the performance of all photochromic lenses are influenced by temperature.


To overcome this challenge, Transitions® XTRActive® lenses use a different formulation of photochromic dyes. *Transitions XTRActive* lenses feature a broad spectrum dye specially designed to react to both UV and visible light for extra darkness outdoors, even in the car. These special molecules collect "extra" energy enabling the lens to darken more, even in the hottest climates.

The logo for Transitions Life360 is centered at the top. It features the word "Transitions" in a blue sans-serif font, followed by a blue circular icon with a white swirl inside. Below this, the word "Life360" is written in a large, black, cursive-style font, with a small "TM" trademark symbol to its upper right. The entire logo is enclosed within a thin black circular border.

Transitions Optical develops photochromic lenses using a proprietary measurement methodology designed to look at the experience of wearers in real world conditions, real locations and real situations. Life360™ includes three types of measurements:



Finally, Transitions Optical has evolved the way photochromics are tested to capture **Real World Measurements**. Real world measurements means exactly that. Products are taken out into the world and the performance is quantified in the many hundreds of different real life conditions. Representing more than a thousand scenarios - combining temperatures, angles of light, UV, weather conditions, time of day, time of year and geographics to reflect how products perform in all situations.




- The visible light range is the portion that allows the eye to see color.
- The UV region, particularly the UVA portion, is the energy range that causes the chemical reaction in photochromic lenses.



SHIFTING INTO THE VISIBLE REGION

If we look at the absorbance* of the dyes used in *Transitions* lenses, these dyes in their colorless form cease to absorb beyond ~410nm. In contrast, the dyes used in *Transitions XTRActive* lenses absorb visible light, out to ~430nm.



The outdoor solar irradiance spectrum has an approximately two-fold increase in intensity at ~400nm vs. ~390nm. Because the absorbance spectra of the *Transitions XTRActive* lens permits this energy to be absorbed and the photochromic rearrangement reaction to occur, the lens is able to darken more as it captures a larger quanta of energy when compared with the *Transitions Signature VII* lens that has less absorbance in this same region of the EM spectrum.

THREE WAYS THIS UNIQUE FORMULA DELIVERS BENEFIT TO YOUR PATIENTS:

1. Designed for extra protection from light, even in the brightest sun and hottest conditions.
2. A hint of tint helps protect eyes from harsh indoor lighting, helping to reduce eye strain and fatigue.
3. Achieves up to category 2 darkness behind the windshield to protect eyes from sunlight while driving.¹

*See glossary.




MORE ACTIVATION: DARKER FOR EXTRA LIGHT PROTECTION OUTDOORS

As with all photochromic technology, *Transitions XTRActive* lenses are influenced by temperature. This is because heat reforms the chemical bond of photochromic molecules - hot temperatures drive the lenses back to clear; and the colder the temperature, the darker the lenses will get. This temperature factor has been minimized in *Transitions XTRActive* lenses.

The visible light activation allows *Transitions XTRActive* lenses to become even darker in both average temperatures and in hot temperatures than traditional photochromic lenses.

Life360 data shows that *Transitions XTRActive* lenses in both gray and brown are darker than *Transitions Signature VII* lenses. Live wearer testing shows that 4 out of 5 wearers are satisfied with the level of darkness outdoors. ²

DARKER IN HOT CLIMATES¹
AVERAGE DARKNESS IN 90°F/32°C



**4 OUT OF 5 WEARERS ARE SATISFIED WITH
THE LEVEL OF DARKNESS OUTDOORS.**



**3 OUT OF 4 WEARERS
ARE SATISFIED WITH
INDOOR CLARITY**

06

A circular graphic with a white background. A thick black horizontal line is at the top. In the center, the text "ACTIVATION BEHIND THE WINDSHIELD" is written in a bold, black, sans-serif font. A vibrant rainbow beam of light enters from the bottom left, angled upwards towards the center. The beam is composed of multiple parallel lines of color, transitioning from purple on the left to red on the right. The text "UV RAYS" is written in a small, black, sans-serif font at the bottom left of the rainbow beam. The entire graphic is framed by a thin black circular border. On the left side, there is a small white circle with four thin black lines radiating from it. On the right side, there are several thin black lines radiating from the outer edge of the circle.

ACTIVATION BEHIND
THE WINDSHIELD

UV RAYS

Additionally, the roof of a vehicle substantially reduces the available irradiance intensity reaching the driver opposed to the intensity present outside in full sunlight.

For these reasons, *Transitions Signature VII* lenses don't offer significant darkness in the car. However, the absorbance spectrum (or, the visible light activation) of the *Transitions XTRActive* lens permits a larger quantum of energy to be absorbed and makes it possible for the lens to darken and activate behind the windshield of a car during the day. In car activation is dependent on temperature, driver position, windshield and lateral window transmission properties.





-

In spite of these excellent wearer results, *Transitions XTRActive* lenses are not a replacement for sunglasses - especially polarized sunglasses which block blinding glare (bright sunlight reflected off surfaces such as the road). They will, however, benefit patients more than a clear lens would in the car.




- even though *Transitions XTRActive* lenses don't block blinding glare behind the windshield - they would still benefit patients more than a clear lens would in the car.

AVERAGE OF REAL WORLD MEASUREMENTS WHEN DRIVING¹

YOUR OPPORTUNITY

Transitions XTRActive lenses offer patients extra protection from light indoors, outdoors and even in the car. Anywhere they go, *Transitions XTRActive* lenses are always working to protect your patients' eyes from the brightest sun and harsh artificial light. *Transitions XTRActive* lenses provide unique benefits that can help you satisfy more of your patients' needs.



**3 OUT OF 4 EYECARE
PROFESSIONALS AGREE³**

Transitions XTRActive lenses are a great option for patients new to the photochromic category.



* GLOSSARY

(IN ORDER OF APPEARANCE)

Irradiance: Irradiance - or light energy from the sun - is the power of electromagnetic radiation per unit area incident on a surface. It is measured in W / m².

Absorbance: The absorbance (also called optical density) of a material is a logarithmic ratio of the radiation falling upon a material, to the radiation transmitted through a material. It is unitless.

SOURCES


¹ Transitions Optical, Inc. Real World Measurements, (2012 - 2013)

² Double blind wearers tests conducted in Brazil (Expertise, 2014), Spain (Ifop, 2014), and US (Pinnacle Marketing, 2009)

³ Transitions XTRActive lenses Eyecare Professional Claims Study conducted in France (MSW-ARS, 2014)

EQUIVALENCE TABLE BETWEEN
ABSORBANCE AND TRANSMISSION

ABSORBANCE	TRANSMISSION
0	100%
0.1	79%
0.25	56%
0.5	32%
0.75	18%
1	10%
2	1%
3	0.1%



VISIT US FOR MORE INFORMATION

 [FACEBOOK.COM/TRANSITIONSOPTICAL](https://www.facebook.com/TRANSITIONSOPTICAL)



**Like us for exclusive
Tech Talk content.**

www.Transitions.com

Transitions, the swirl and XTRActive are registered trademarks and *Transitions Signature* and *Life360* are trademarks of Transitions Optical, Inc. ©2015 Transitions Optical, Inc.
Photochromic performance is influenced by temperature, UV exposure and lens material. LTRN200506 SHK/TOI 4/15